

VarioCup, Cemented




VarioCup, Cemented

highly crosslinked polyethylene to reduce wear rate



Stainless steel ring works as x-ray marker



Leverage traditional Muller Cup design
Neutral and hooded versions
Wide size range - 47-59mm
3mm cement mantle thickness, 1.5mm cement thickness at rim region.(the labelled diameter includes 3mm cement mantle thickness)

Same Tray Instrument Sets for VarioCup, only additional cup holder is needed.

Cup Head Compatibility Chart

Cup Size in mm	47	49	51	53	55	57	59
Head Size	28mm						
			32mm				

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.